

All India Institute of Medical Sciences
Saket Nagar Bhopal

Pre-Bid Clarification and Corrigendum

Dated : 03 Feb. 2017

With reference to E Tender No. 004 AIIMS BHOPAL/MC/CS -ANATOMY/HISTOLOGY LAB EQUIPMENT/2016-17 /01 Dated :16.01.2017 for Notice .Inviting Tender for Supply. Installation and support of Anatomy/Histology Lab. Equipments.

Subsequent to Pre-Bid queries sought for clarification in Pre-Bid meeting or queries received through e-mail/hard copy ,the following clarification to Pre-Bid queries are circulated for the wider information to the prospective tenders and requested them to download Pre-Bid clarification and submit the same along with their "Technical Bid" envelope .

Queries raised by

1.M/S ESCO Biotech Pvt. Ltd., Mumbai

S.No	Tender Specification	Query	Clarification
1	Item no. 7 CO2 Incubator Low Water Alarm Indication	Low water indication comes in water jacked CO2 Incubators. In CO2 incubator, RH level to be monitored so equipment should have RH display	Low water alarm indication/RH display/Humidity indicator
2	Item no. 7 CO2 Incubator Communication protocols HL-7 for networked environments to HIS	This is not needed for CO2 . HL 7 allows the integration of various applications, such as bedside terminals, radiological imaging station, hospital census, order entries and patient accounting, into one system . This is useful for diagnostic equipment for pathology.	Is retained
3	Item No. 9 Biological Safety Cabinet Dimensions (cabinet Size) 1200 to 1800 feet	This should be in mm	<ul style="list-style-type: none">• Is typographical error• Dimensions (cabinet Size) 1200 to 1800 mm
4	Item no. 12 Deep Freezer – 80 c Access card option for multiple users is highly desirable	Only needed for pharmaceutical companies where only one or two people uses the freezer, Most of the freezers comes with password protection to protect from unauthorized access.	Access protection by password/Access card/Front lock & key

Corrigendum

S.No.	Specification	Recommendation
1	CO2 Incubators (Hypoxia & Normoxias)	<ul style="list-style-type: none">• Remove Hypoxia & Normoxias• Typographical error• Retain only CO2 Incubator

S. K. Singh
3/2/2017
for
Head
Deptt. Of Anatomy
AIIMS Bhopal